

SOX11 Rabbit mAb

货号: B30057

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	46kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,SH-SY5Y,Neuro-2a,Y79,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SOX11.
序列	MVQQAESLEAESNLPREALDTEEGEFMACSPVALDESDPDWCKTASGHIKRPMAFMVWSKIERRKIMEQSPDMHNAEISKRLGKRWKMLKDSEKIPFIR

靶点信息

研究背景	This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The protein may function in the developing nervous system and play a role in tumorigenesis.
基因ID	6664

基因名	SOX11
Swiss	P35716
别名	SOX11;MRD27

产品验证



Western blot analysis of SOX11 expressed in Mouse brain,Rat brain,BT-474 using SOX11 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn