

# Cdc25C Rabbit mAb

货号: B30066

产品信息

|       |                                                                                                                                          |
|-------|------------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                                                    |
| 宿主    | Rabbit                                                                                                                                   |
| 克隆性   | Monoclonal                                                                                                                               |
| 预测反应  |                                                                                                                                          |
| 应用    | WB IHC IF/ICC IP FC                                                                                                                      |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 45kDa/48kDa/50kDa/53kDa                                                                                                                  |
| 实测分子量 | 60kDa                                                                                                                                    |
| 形式    | Liquid                                                                                                                                   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                                            |
| 偶联物   | Unconjugated                                                                                                                             |
| 阳性对照  | Mouse lung,Mouse thymus,Rat thymus                                                                                                       |
| 细胞定位  | Nucleus                                                                                                                                  |
| 纯化    | Affinity purification                                                                                                                    |

抗原信息


|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human Cdc25C.                                                          |
| 序列   | SFRSNQRKMLNLLLERDTSFTVCPDVPRTVPVGKFLGDSANLSILSGGTPKRCLDLNLSSGEITATQLTTSADLDETGHLD<br>SSGLQEVHLAGMNHQHLMKCSPAQLLCST |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                             |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a conserved protein that plays a key role in the regulation of cell division. The encode d protein directs dephosphorylation of cyclin B-bound CDC2 and triggers entry into mitosis. It also suppre sses p53-induced growth arrest. Multiple alternatively spliced transcript variants of this gene have been d escribed. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                      |
|-------|----------------------|
| 基因ID  | 995                  |
| 基因名   | CDC25C               |
| Swiss | P30307               |
| 别名    | CDC25C;CDC25;PPP1R60 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)