

Somatostatin Receptor 2 Rabbit mAb

货号: B30088

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	40kDa/41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,C6,SH-SY5Y,Neuro-2a,MDA-MB-231,Mouse thymus
细胞定位	Cell membrane,Cytoplasm,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Somatostatin Receptor 2.
序列	YANSCANPILYAFLSDNFKSFQNVLCCLKVSGTDDGERSDSKQDKSRLNETTETQRTLLNGDLQTSI

靶点信息

研究背景	Somatostatin acts at many sites to inhibit the release of many hormones and other secretory proteins. The biologic effects of somatostatin are probably mediated by a family of G protein-coupled receptors that are expressed in a tissue-specific manner. SSTR2 is a member of the superfamily of receptors having seven transmembrane segments and is expressed in highest levels in cerebrum and kidney.
基因ID	6752

基因名	SSTR2
Swiss	P30874
别名	SSTR2;somatostatin receptor 2;SS 2 R;SS2R;SRIF 1;SS2 R;SRIF-1

产品验证



Western blot analysis of Somatostatin Receptor 2 expressed in Mouse thymus,MCF7,C6 using Somatostatin Receptor 2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn