

CETP Rabbit mAb

货号: B30095

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	47kDa/54kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HepG2,A375,Mouse liver,Rat liver
细胞定位	Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CETP.
序列	LMSLMGDEFKAVLETWGFNTNQEIFQEVVGGFPSQAQVTVHCLKMPKISCQNKGVVVNSSVMVKFLFPRPDQQHSVA YTFEEDIVTTVQASYSKKFLSLDFQITPKTVSNLTSSSESQSFLQSMITAVGIPEVMSRLEVVFTALMNSKGVSLFDII NPEIITRDGFLLLQMDFGFPEHLLVDFLQSL

靶点信息

研究背景	The protein encoded by this gene is found in plasma, where it is involved in the transfer of cholesteryl ester from high density lipoprotein (HDL) to other lipoproteins. Defects in this gene are a cause of hyperalphalipoproteinemia 1 (HALP1). Two transcript variants encoding different isoforms have been found for this gene.
基因ID	1071

基因名	CETP
Swiss	P11597
别名	CETP;BPIFF;HDLQC10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn