

BMP7 Rabbit mAb

货号: B30103

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	49kDa
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,293T,MCF7
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human BMP7.
序列	PHLQGKHNSAPMFMLDLYNAMAVEEGGGPGGQGFSPYKAVFSTQGPPLASLQDSHFLTDADMVMSFVNLVEHDKEFF HPRYHHREFRFDLSKIPEGEAVTAAEFRIYKDYIRERFDNETFRISVYQVLQEHLGRESDLFLLDSTRTLWASEEGWLVDITA TSNHVVVNPRHNLGLQLSVET

靶点信息

研究背景	This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMA D family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which plays a role in bone, kidney and brown adipose tissue development. Additionally, this protein induces ectopic bone formation and may promote fracture healing in human patients.
------	--

基因ID	655
基因名	BMP7
Swiss	P18075
别名	BMP7;OP-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn