

GATA3 Rabbit mAb

货号: B30108

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC CHIP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 CHIP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	47kDa/48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human GATA3.
序列	MEVTADQPRWVSHHHPAVLNGQHPDTHHPGLSHSYMDAAQYPLPEEVDVLFNIDGQGNHVPPYYGNSVRATVQRYPPT HHGSQVCRPPLLHGSLPWLDGGKALGSHHTASPWNLSPFSKTSIIHHGSPGPLSVYPPASSSSLSGGHASPHLFTFPPTPP KDVSPDPSLSTPGSAGSARQDEKECLKYQVPLPDSMKLESSHSRGSM TALGGASSSTHHPITTYPPYVPEYSSGLFPPSSL LGGSP TGF GCKSRPKARSSTG

靶点信息

研究背景	This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.
基因ID	2625
基因名	GATA3
Swiss	P23771
别名	GATA3;HDR;HDRS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn