

PFDN5 Rabbit mAb

货号: B30125

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	7kDa/12kDa/17kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,U-87MG,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PFDN5.
序列	MAQSINITELNLPQLEMLKNQLDQEVEFLSTSIAQLKVVTQKYVEAKDCLNVLNKSNEGKELLVPLTSSMYVPGKLHDVEHV LIDVGTGYVVEKTAEDAKDFFKRKIDFLTKQMEKIQPALQEKHAMKQAVMEMMSQKIQQLTALGAAQATAKA

靶点信息

研究背景	This gene encodes a member of the prefoldin alpha subunit family. The encoded protein is one of six sub units of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptid es, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils. The encoded protein may also repre ss the transcriptional activity of the proto-oncogene c-Myc. Alternatively spliced transcript variants encodi ng different isoforms have been described.
------	--

基因ID	5204
基因名	PFDN5
Swiss	Q99471
别名	PFDN5;MM-1;MM1;PFD5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn