

NEFH Rabbit mAb

货号: B30163

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	200kDa
实测分子量	200kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse brain,Rat brain
细胞定位	axon,cytoplasm,cytoskeleton,neurofilament,postsynaptic density,postsynaptic intermediate filament cyto skeleton,Schaffer collateral - CA1 synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NEFH.
序列	AKSPAKEEAKSPPEAKSPEKEEAKSPA EVKSPEKAKSPAKEEAKSPA EAKSPEKAKSPVKEEAKSPA EAKSPVKEEAKSPA E VKSPEKAKSPTKEEAKSPE

靶点信息

研究背景	Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and functionally maintain neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene encodes the heavy neurofilament protein. This protein is commonly used as a biomarker of neuronal damage and susceptibility to amyotrophic lateral sclerosis (ALS) has been associated with mutations in this gene. [provided by RefSeq, Oct 2008]
------	--

基因ID	4744
基因名	NEFH
Swiss	P12036
别名	CMT2CC;NFH;NEFH

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn