

# UBE2B Rabbit mAb

货号: B30166

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	17kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,SW480,Jurkat,NIH/3T3,Mouse heart,Mouse brain,Mouse skeletal muscle,Rat heart,Rat brain
细胞定位	Cell membrane,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human UBE2B.
序列	MSTPARRRLMRDFKRLQEDPPVGVSGAPSENNIMQWNAVIFGPEGTPFEDGTFKLVI EFSEEYPNKPPTVRFLSKMFHPNV YADGSICLDILQNRWSPTYDVSSILTSIQSLLDEPNPNSPAN SQAAQLYQENKREYEKRVS AIVEQSWNDS

## 靶点信息

研究背景	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or s hort-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-ac tivating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. Thi s gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is required for pos t-replicative DNA damage repair. Its protein sequence is 100% identical to the mouse, rat, and rabbit ho mologs, which indicates that this enzyme is highly conserved in eukaryotic evolution.
------	---

基因ID	7320
基因名	UBE2B
Swiss	P63146
别名	UBE2B;E2-17kDa;HHR6B;HR6B;RAD6B;UBC2

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)