

Phospho-JunD (S255) Rabbit mAb

货号: B30185

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	35kDa
实测分子量	38,42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,293T,HepG2,Mouse kidney,Mouse liver,Mouse pancreas,Mouse testis,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Phospho-JunD (S255).
序列	METPFYGDEALSGLGGGASGSGGSFASPGRLFPGAPPTAAAGSMMKKDALTLSLSEQVAAALKPAAAPPPTPLRADGAP SAAPPDGLLASPDGLLLKLASPELERLIQSNGLVTTTPTSSQFLYPKVAASEEQEFAEGFVKALEDLHKQNQLGAGAAAAA AAAAAGGPSGTATGSAPPGELAPAAAAPE

靶点信息

研究背景	The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216).
基因ID	3727

基因名	JUND
Swiss	P17535
别名	JUND;AP-1

产品验证



Western blot analysis of Phospho-JunD (S255) expressed in Mouse kidney,Rat liver,HeLa using Phospho-JunD (S255) Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn