

Tetranectin Rabbit mAb

货号: B30208

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	22kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Tetranectin.
序列	EPPTQKPKKIVNAKKDVVNTKMFEELKSRLDTLAQEVALLEQQALQTVCLKGTVHMKCFLAFTQTKTFHEASEDCISRG GTLGTPQTGSENDALYEYLRQSVGNEAEIWLGLNDMAAEGTWVDMTGARIAYKNWETEITAQPDGGKTENCAVLSGAA NGKWFDKRCRDQLPYICQFGIV

靶点信息

研究背景	Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecule s destined for exocytosis.
基因ID	7123

基因名	CLEC3B
Swiss	P05452
别名	CLEC3B;TN;TNA

产品验证



Western blot analysis of Tetranectin expressed in HUVEC,mouse liver,rat heart using Tetranectin Rab bit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure ti me: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn