

CPS1 Rabbit mAb

货号: B30209

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	116kDa/164kDa/165kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U2OS
细胞定位	Mitochondrion,Nucleus,nucleolus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human CPS1.
序列	IKNVVSGKTSACFEPSLDYMVTKIPRWDLDRFHGTSSRIGSSMKSVGEVMAIGRTFEESFQKALRMCHPSIEGFTPRLPMN KEWPSNLDLRKELSEPSSTRIYAIKAIDDNMSLDEIEKLTIDKWFLYKMRDILNMEKTLKGLNSESMTTEETLKRAKEIGFS DKQISKCLGLTEAQTRELRLKKNHPWVKQIDTLAAEYPSVTNYLYVTYNGQEHDVN

靶点信息

研究背景	The mitochondrial enzyme encoded by this gene catalyzes synthesis of carbamoyl phosphate from ammonia and bicarbonate. This reaction is the first committed step of the urea cycle, which is important in the removal of excess urea from cells. The encoded protein may also represent a core mitochondrial nucleoid protein. Three transcript variants encoding different isoforms have been found for this gene. The shortest isoform may not be localized to the mitochondrion. Mutations in this gene have been associated with carbamoyl phosphate synthetase deficiency, susceptibility to persistent pulmonary hypertension, and susceptibility to venoocclusive disease after bone marrow transplantation.
基因ID	1373
基因名	CPS1
Swiss	P31327
别名	CPS1;CPSASE1;PHN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn