

MUC5AC Rabbit mAb

货号: B30214

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	585kDa
实测分子量	220-600kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse stomach
细胞定位	cytoplasm,extracellular exosome,extracellular space,Golgi lumen,mucus layer,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human MUC5AC.
序列	PLSCSLDEARMSKDGCCRFPPPPPPYQNQSTCAVYHRSIIQQGCSSEPVRLAYCRGNCGDSSSMYSLEGNTVEHRC QCCQELRTSLRNVTLHCTDGSSRAFSYTEVEECGCMGRRCPAPGDTQHSEEAPEPSQEAESGSWERGVVPSPMH

靶点信息

研究背景	Gel-forming glycoprotein of gastric and respiratory tract epithelia that protects the mucosa from infection and chemical damage by binding to inhaled microorganisms and particles that are subsequently removed by the mucociliary system. Interacts with H.pylori in the gastric epithelium, Barrett's esophagus as well as in gastric metaplasia of the duodenum (GMD).
基因ID	4586

基因名	MUC5AC
Swiss	P98088
别名	MUC5AC;MUC5;TBM;leB;mucin;mucin-5AC;MUC5A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn