

HOXA13 Rabbit mAb

货号: B30220

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	39kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,HT-29
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human HOXA13.
序列	GPAGPAGAEAAKQCSPCSAAAQSSSGPAALPYGYFGSGYPCARMGPHNPAIKSCAQPASAAAAAFADKYMDTAGPA AEEFSSRAKEFAFYHQGYAAGPYHHHQMPGYLDMPVVPGLGGPGESRHEPLGLPMESYQPWALPNGWNGQMYCPKE QAQPPHLWKSTLPDVVSHPSDASSYR

靶点信息

研究背景	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Expansion of a polyalanine tract in the encoded protein can cause hand-foot-uterus syndrome, also known as hand-foot-genital syndrome.
------	--

基因ID	3209
基因名	HOXA13
Swiss	P31271
别名	HOXA13;HOX1;HOX1J

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn