

Hemopexin Rabbit mAb

货号: B30223

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	51kDa
实测分子量	70-75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,HepG2,BT-474,Mouse liver,Mouse pancreas
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Hemopexin.
序列	TPLPPTSAHGNVAEGETKPDPDVTERCSDGWSFDATTLDDNGTMLFFKGEFVWKSHKWDRELISERWKNFPSPVDAAF RQGHNSVFLIKGDKVWVYPPEKKEKGYPKLLQDEFPGIPSPLDAAVECHRGECQAEGVLFFQGDREWFWDLATGTMKER SWPAVGNCSSALRWLGRIYCFQGNQFLRFDPPVRGEVPPRYPRDVRDYFMPCPGRGHGHRNGTGHG

靶点信息

研究背景	This gene encodes a plasma glycoprotein that binds heme with high affinity. The encoded protein is an acute phase protein that transports heme from the plasma to the liver and may be involved in protecting cells from oxidative stress.
------	--

基因ID	3263
基因名	HPX
Swiss	P02790
别名	HPX;HX;hemopexin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn