

Caldesmon Rabbit mAb

货号: B30231

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	61kDa/62kDa/64kDa/65kDa/93kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431
细胞定位	Cytoplasm,cytoskeleton,myofibril
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Caldesmon.
序列	MDDFERRRELRRQKREEMRLEAERIAYQRNDDDEEEAARERRRRRARQERLRQKQEEESLGQVTDQVEVNAQNSVPDEEA KTTTTNTQVEGDDEAAFLERLARREERRQKRLQEALERQKEFDPTITDASLSLPSRRMQNDTAENETTEKEEKSESQRERY EIEETETVTKSYQKNDWRDAEENKKEDKEKEEEEEEEKPKRGSIGENQIKDEKIKDKKEPKKEEVKSFMDRKKGFTEVKSQN

靶点信息

研究背景	This gene encodes a calmodulin- and actin-binding protein that plays an essential role in the regulation of smooth muscle and nonmuscle contraction. The conserved domain of this protein possesses the binding activities to Ca(2+)-calmodulin, actin, tropomyosin, myosin, and phospholipids. This protein is a potent inhibitor of the actin-tropomyosin activated myosin MgATPase, and serves as a mediating factor for Ca(2+)-dependent inhibition of smooth muscle contraction. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms.
基因ID	800
基因名	CALD1
Swiss	Q05682
别名	CALD1;CDM;H-CAD;HCAD;L-CAD;LCAD;NAG22;caldesmon

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn