

# CDC7 Rabbit mAb

货号: B30275

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB FC                                                                                         |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>FC:</b> 1:20 - 1:50                                           |
| 理论分子量 | 63kDa                                                                                         |
| 实测分子量 | 70kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | SH-SY5Y,HepG2,THP-1                                                                           |
| 细胞定位  | Nucleus                                                                                       |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human CDC7 .                                                                                                                                                                                                                                                                |
| 序列   | EGSVGLSVQRSVFGERNFNHSSISHESPAVKLMKQSKTVDVLSRKLATKKKAISTKVMNSAVMRKTASSCPASLTCDCYA<br>TDKVCISICLSRRQQVAPRAGTPGFRAPEVLTCPNQTTAIDMWSAGVIFLSLLSGRYPFYKASDDLTALAQIMTIRGSRETI<br>QAAKTFGKSILCSKEVPAQDLRKLCELRGMDSSTPKLTSDIQGHASHQPAISEKTDHKASCLVQTTPPGQYSGNSFKKGD<br>SNSCEHCFDEYNTNLEGWNEVPDEAYDLLDKLLDLPASRITAEAEALLHPFFKDMSL |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a cell division cycle protein with kinase activity that is critical for the G1/S transition. T he yeast homolog is also essential for initiation of DNA replication as cell division occurs. Overexpression of this gene product may be associated with neoplastic transformation for some tumors. Multiple alternati vely spliced transcript variants that encode the same protein have been detected. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                |
|-------|--------------------------------|
| 基因ID  | 8317                           |
| 基因名   | CDC7                           |
| Swiss | O00311                         |
| 别名    | CDC7;CDC7L1;HsCDC7;Hsk1;huCDC7 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)