

FHL1 Rabbit mAb

货号: B30350

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	22kDa/31kDa/33kDa/34kDa/36kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,22Rv1,MCF7,Mouse kidney
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human FHL1.
序列	RDPLQGKKYVQKDGHHCLKCFDKFCANTCVEC RKPIGADSKVEHYKNRFWHDTCFRCAKCLHPLANETFVAKDNKILC NKCTTREDSPKCKGCFKAIVAGDQNV EYKGTVWHKDCFTCSNCKQVIGTGSFFPKGEDFYCVTCHETKFAKHCVKCNKAI TSGGITYQDQPWHADCFVCVTC SKKLAGQRFTAVEDQYYCVD CYKNFVAKKCAGCKNPITGFGKGSSV VAYEGQSWHD YCFHCKKCSVNLANKRFVFHQEQVYCPDCAKKL

靶点信息

研究背景	This gene encodes a member of the four-and-a-half-LIM-only protein family. Family members contain two highly conserved, tandemly arranged, zinc finger domains with four highly conserved cysteines binding a zinc atom in each zinc finger. Expression of these family members occurs in a cell- and tissue-specific mode and these proteins are involved in many cellular processes. Mutations in this gene have been found in patients with Emery-Dreifuss muscular dystrophy. Multiple alternately spliced transcript variants which encode different protein isoforms have been described.
基因ID	2273
基因名	FHL1
Swiss	Q13642
别名	FHL1;FHL-1;FHL1A;FHL1B;FLH1A;KYOT;RBMX1A;RBMX1B;SLIM;SLIM-1;SLIM1;SLIMMER;XMPMA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn