

LMO4 Rabbit mAb

货号: B30377

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	17kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human LMO4.
序列	MVNP GSSSQPPV TAGSLSWKRCAGCGGKIADRFLLYAMD SYWHSRCLKSCCQAQLGDIGTSCYTKSGMILCRNDYIR LFGNSGAC SACGQSIPASELVMRAQGNVYHLKCFT CSTCRNRLVPGDRFHYINGSLFCEHDRPTALINGHLNSLQSNPLLP DQKVC

靶点信息

研究背景	This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeo domain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.
基因ID	8543

基因名	LMO4
Swiss	P61968
别名	LMO4

产品验证



Western blot analysis of LMO4 expressed in human breast cancer,mouse brain using LMO4 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn