

JunD Rabbit mAb

货号: B30395

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	35kDa
实测分子量	38,42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,293T,HepG2,Mouse kidney,Mouse liver,Mouse pancreas,Mouse testis,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human JunD.
序列	METPFYGD E ALSGLG G ASGSGGSFASPGRLFPGAPPTAAAGSMMKKDAL T LSLSEQVAAALKPAAAPPPTPLRADGAPSAAPPDGLLASPD L GLLKLASPELERLI I QSNGLVTTTPTSSQFLYPKVAASEEQEFAEGFVKAL E DLHKQNQLGAGAAAAA AAAAAGGPGSGTATGSAPP G ELAPAAAAPE

靶点信息

研究背景	The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216).
------	--

基因ID	3727
基因名	JUND
Swiss	P17535
别名	JUND;AP-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn