

HSPA4 Rabbit mAb

货号: B30412

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	110kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,NIH/3T3,K-562,C6,Mouse testis,Mouse brain,Rat brain
细胞定位	cytosol,extracellular exosome,mitochondrion,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human HSPA4.
序列	EELGKQIQQYMKIISSEFNKEDQYDHLDAADMTKVEKSTNEAMEWMNNKLNQNKQSLTMDPVVKSKEIEAKIKELTSTC SPIISKPKPKVEPPKEEQKNA

靶点信息

研究背景	Predicted to enable ATP binding activity. Involved in chaperone-mediated protein complex assembly and protein insertion into mitochondrial outer membrane. Located in cytosol and extracellular exosome. Implicated in Chagas disease. Biomarker of chronic obstructive pulmonary disease; rheumatoid arthritis; type 2 diabetes mellitus; and ulcerative colitis.
------	--

基因ID	3308
基因名	HSPA4
Swiss	P34932
别名	APG-2; HEL-S-5a; HS24/P52; HSPH2; RY; hsp70; hsp70RY

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn