

CRMP1 Rabbit mAb

货号: B30413

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	62kDa/74kDa
实测分子量	62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat brain
细胞定位	Cytoplasm,centrosome,cytoskeleton,microtubule organizing center,spindle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CRMP1.
序列	NVNKGMGRFIPRKAFPEHLYQRVKIRNKVFGGLQGVSRGMYDGPVYEVPA TP KYATPAPSAKSSPSKHQPPPIRNLHQSNF S

靶点信息

研究背景	This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants.
------	---

基因ID	1400
基因名	CRMP1
Swiss	Q14194
别名	CRMP1;CRMP-1;DPYSL1;DRP-1;DRP1;ULIP-3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn