

beta 2 Microglobulin Rabbit mAb

货号: B30417

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC IP
推荐浓度	WB: 1:500 - 1:2000 ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	13kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,COS-7,HeLa,K-562,Mouse lung,Mouse liver,Mouse spleen
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human beta 2 Microglobulin.
序列	IQRTPKIQVYSRHPAENGKSNFLNCYVSGFHPSDIEVDLLKNGERIEKVEHSDLSFSKDWSFYLLYYTEFTPTEKDEYACRV NHVTLSQPKIVKWDRDM

靶点信息

研究背景	This gene encodes a serum protein found in association with the major histocompatibility complex (MHC) class I heavy chain on the surface of nearly all nucleated cells. The protein has a predominantly beta-plea ted sheet structure that can form amyloid fibrils in some pathological conditions. The encoded antimicrob ial protein displays antibacterial activity in amniotic fluid. A mutation in this gene has been shown to resu lt in hypercatabolic hypoproteinemia.
------	--

基因ID	567
基因名	B2M
Swiss	P61769
别名	B2M;IMD43

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn