

WNT2 Rabbit mAb

货号: B30420

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,NIH/3T3,Rat brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human WNT2.
序列	SWWYMRATGGSSRVMCDNVPGLVSSQRQLCHRHPDVMRAISQGVAEWTAECQHQFRQHRWNCNTLDRDHSLFGRV LLRSSRESAFVYAISSAGVVFAITRACSQGEVKSCSDPKKMGSAKDSKGIFDWGGCSDNIDYGIKFARAFVDAKERKKG DARALMNLHNNRAGRKAVKRFLKQECKCHGVSGSCTLRTCWLAMADFRKTGDYLRKYNGAIQVVMNQDGTGFTVAN ERFKKPTKNDLVYFENSPDYCIRDREAGSLGTAGRVCNLTSRGMDSCEVMCCGRGYDTSHVTRMTKCGCKFWCCAVR CQDCLEALDVHTCKAPKNADWTTAT

靶点信息

研究背景	This gene is a member of the WNT gene family. The WNT gene family consists of structurally related gene s which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in s everal developmental processes, including regulation of cell fate and patterning during embryogenesis. Al ternatively spliced transcript variants have been identified for this gene.
------	---

基因ID	7472
基因名	WNT2
Swiss	P09544
别名	WNT2;INT1L1;IRP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn