

# PBRM1 Rabbit mAb

货号: B30429

## 产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Monoclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC FC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 193kDa                                                                                                         |
| 实测分子量 | 193kDa                                                                                                         |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                  |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | HeLa,NIH/3T3,Jurkat,Rat brain                                                                                  |
| 细胞定位  | Nucleus                                                                                                        |
| 纯化    | Affinity purification                                                                                          |

## 抗原信息


|      |                                                                                                                                                           |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human PBRM1.                                                                                                  |
| 序列   | VLEAREPGSGRRLCDLFMVKPSKKDYPDYKIIILEPMDLKIIEHNIRNDKYAGEEGMIEDMKLMFRNARHYNEEGSQVYND<br>AHILEKLLKEKRKELGPLDDDDMASPKLKLRSKGISPKSKYMTMPQQKLNEVYEAVKNYTDKRGRRRL |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                               |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This locus encodes a subunit of ATP-dependent chromatin-remodeling complexes. The encoded protein has been identified as an integral component of complexes necessary for ligand-dependent transcriptional activation by nuclear hormone receptors. Mutations at this locus have been associated with primary clear cell renal cell carcinoma. [provided by RefSeq, Feb 2012] |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |             |
|-------|-------------|
| 基因ID  | 55193       |
| 基因名   | PBRM1       |
| Swiss | Q86U86      |
| 别名    | BAF180; PB1 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)