

# MFF Rabbit mAb

货号: B30449

## 产品信息

|       |                                                                                                                                     |
|-------|-------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                     |
| 宿主    | Rabbit                                                                                                                              |
| 克隆性   | Monoclonal                                                                                                                          |
| 预测反应  |                                                                                                                                     |
| 应用    | WB                                                                                                                                  |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                                                           |
| 理论分子量 | 34kDa                                                                                                                               |
| 实测分子量 | 34kDa                                                                                                                               |
| 形式    | Liquid                                                                                                                              |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                                       |
| 偶联物   | Unconjugated                                                                                                                        |
| 阳性对照  | C2C12,C6,Rat brain                                                                                                                  |
| 细胞定位  | Cytoplasmic vesicle,Mitochondrion outer membrane,Peroxisome,Single-pass type IV membrane protein,secretory vesicle,synaptic vesicle |
| 纯化    | Affinity purification                                                                                                               |

## 抗原信息


|      |                                                                                                      |
|------|------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human MFF.                                               |
| 序列   | TPSPQQARVCPHMLPEDGANLSSARGILSLIQSSTRRAYQQILDVLDENRRPVLRGGSAAATSNPHHDNVRYGISNIDTTIEGTSDDLTVVDAASLRRQI |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This is a nuclear gene encoding a protein that functions in mitochondrial and peroxisomal fission. The encoded protein recruits dynamin-1-like protein (DNM1L) to mitochondria. There are multiple pseudogenes for this gene on chromosomes 1, 5, and X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013] |
| 基因ID | 56947                                                                                                                                                                                                                                                                                                                                                 |

|       |                       |
|-------|-----------------------|
| 基因名   | MFF                   |
| Swiss | Q9GZY8                |
| 别名    | C2orf33; EMPF2; GL004 |

产品验证



Western blot analysis of MFF expressed in C2C12,C6,Rat brain using MFF Rabbit mAb at 1:1000. Second ary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buff er: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)