

RAB11A Rabbit mAb

货号: B30489

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	24kDa
实测分子量	24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,SH-SY5Y,Mouse testis,Mouse brain,Rat testis,Rat brain,Rat spleen
细胞定位	Cell membrane,Cleavage furrow,Cytoplasmic side,Cytoplasmic vesicle,Lipid-anchor,Peripheral membrane protein,Recycling endosome membrane,phagosome,phagosome membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RAB11A.
序列	ENVERWLKELRDHADSNIIVMLVGNKSDLRHLRAVPTDEARAFAEKNGLSFIETSALDSTNVEAAFQILTEIYRIVSQKQMS DRRENDMSPSNNVPIHV

靶点信息

研究背景	The protein encoded by this gene belongs to the Rab family of the small GTPase superfamily. It is associated with both constitutive and regulated secretory pathways, and may be involved in protein transport. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2011]
基因ID	8766

基因名	RAB11A
Swiss	P62491
别名	YL8

产品验证



Western blot analysis of RAB11A expressed in Mouse testis,Rat testis,HeLa using RAB11A Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn