

SRP54 Rabbit mAb

货号: B30523

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	50kDa/55kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse brain,Mouse heart,Mouse liver,Mouse lung
细胞定位	Cytoplasm,Nucleus speckle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SRP54.
序列	VIVTKLDGHAKEGGGALSAVAATKSPIIFIGTGEHIDDFEPFKTQPFISKLLGMGDIEGLIDKVNELKLDDNEALIEKLKHGQFT LRDMYEQFQNMKMGPFQILGMIPGFGTDFMSKGNEQESMARLKKLMTIMDSMNDQELDSTDGAKVFSKQPGRIQRV ARGSGVSTRDVQELLTQYTKFAQMVKKMGGIKGLFKGGDMSKNVSQSQMAKLNQQMAKMMDPRVLHHMGGMAGLQ SMMRQFQQGAAGNMKGMMGFNNM

靶点信息

研究背景	Binds to the signal sequence of presecretory protein when they emerge from the ribosomes and transfers them to TRAM (translocating chain-associating membrane protein. Plays a role in proliferation and differe ntiation of granulocytic cells, neutrophils migration capacity and exocrine pancreas development.
------	--

基因ID	6729
基因名	SRP54
Swiss	P61011
别名	SRP54

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn