

Phospho-GSK3 (alpha + beta)(Y216 + Y279) Rabbit mAb

货号: B30532

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	50kDa
实测分子量	51kDa/46kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-431,HeLa,B-cell,Mouse liver
细胞定位	apical dendrite,axon,beta-catenin destruction complex,cytoplasm,cytosol,mitochondrion,nucleus,postsynapse,proximal dendrite
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Phospho-GSK3 (alpha + beta)(Y216 + Y279).
序列	LGTQLPNNRPLPLFNFSAGELSIQPSLNAILIPPHLRSPAGTTTLTPSSQALTETPTSSDWQSTDATPTLTNSS

靶点信息

研究背景	This gene encodes a multifunctional Ser/Thr protein kinase that is implicated in the control of several regulatory proteins including glycogen synthase, and transcription factors, such as JUN. It also plays a role in the WNT and PI3K signaling pathways, as well as regulates the production of beta-amyloid peptides associated with Alzheimer's disease.
------	---

基因ID	2931
基因名	GSK3A
Swiss	P49840
别名	GSK3 alpha;GSK3A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn