

GLP2 Rabbit mAb

货号: B30577

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa
实测分子量	4kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,293T,Rat spleen
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GLP2.
序列	RSLQDTEEKSRFSASQADPLSDPDQMNEDKRHSQGTFTSDYSKYLSRRAQDFVQWLMNTKRNRRNNIAKRHDEFER HAEGTFTSDVSSYLEGQAAKEFIAWLVKGRGRRDFPEEVAIVEELGRRHADGSFSDENMTILDNLAARDFINWLIQTKITD RK

靶点信息

研究背景	The protein encoded by this gene is actually a preproprotein that is cleaved into four distinct mature peptides. One of these, glucagon, is a pancreatic hormone that counteracts the glucose-lowering action of insulin by stimulating glycogenolysis and gluconeogenesis. Glucagon is a ligand for a specific G-protein linked receptor whose signalling pathway controls cell proliferation. Two of the other peptides are secreted from gut endocrine cells and promote nutrient absorption through distinct mechanisms. Finally, the fourth peptide is similar to glicentin, an active enteroglucagon.
------	---

基因ID	2641
基因名	GCG
Swiss	P01275
别名	GCG;GLP-1;GLP1;GLP2;GRPP;glucagon

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn