

DOK1 Rabbit mAb

货号: B30595

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	65kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,BxPC-3
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human DOK1.
序列	EPKDAWWCQARVKEEGYELPYNPATDDYAVPPPRSTKPLLAPKPQGPAFPEPGTATGSGIKSHNSALYSQVQKSGASGS WDCGLSRVGTDKTGVKSEGST

靶点信息

研究背景	The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosin e kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multipr otein signaling complexes. Several transcript variants encoding different isoforms have been found for thi s gene.
基因ID	1796

基因名	DOK1
Swiss	Q99704
别名	DOK1; P62DOK; pp62; docking protein 1

产品验证



Western blot analysis of DOK1 expressed in MCF7,BxPC-3 using DOK1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn