

# Syntaxin 3 Rabbit mAb

货号: **B30625**

产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Monoclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC FC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 33kDa                                                                                                          |
| 实测分子量 | 33kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                  |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | HT-29,293T,PC-12                                                                                               |
| 细胞定位  | Membrane,Single-pass type IV membrane protein                                                                  |
| 纯化    | Affinity purification                                                                                          |

抗原信息


|      |                                                                                                    |
|------|----------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human Syntaxin 3.                                      |
| 序列   | KHIEEDEVSSADLRIRKSQHSVLSRKFEVMTKYNEAQVDFRERSKGRIQRQLEITGKKTDEELEEMLESGNPAIFTSGIIDSQISKQALSEIEGRHKD |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                   |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The gene is a member of the syntaxin family. The encoded protein is targeted to the apical membrane of epithelial cells where it forms clusters and is important in establishing and maintaining polarity necessary for protein trafficking involving vesicle fusion and exocytosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010] |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |        |
|-------|--------|
| 基因ID  | 6809   |
| 基因名   | STX3   |
| Swiss | Q13277 |
| 别名    | STX3A  |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)