

# NAT10 Rabbit mAb

货号: B30636

产品信息

|       |                                                                                                                                          |
|-------|------------------------------------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                                          |
| 宿主    | Rabbit                                                                                                                                   |
| 克隆性   | Monoclonal                                                                                                                               |
| 预测反应  |                                                                                                                                          |
| 应用    | WB IHC IF/ICC IP FC                                                                                                                      |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 116kDa                                                                                                                                   |
| 实测分子量 | 116kDa                                                                                                                                   |
| 形式    | Liquid                                                                                                                                   |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                                            |
| 偶联物   | Unconjugated                                                                                                                             |
| 阳性对照  | HeLa,K-562,SH-SY5Y,Mouse testis                                                                                                          |
| 细胞定位  | Nucleus,nucleolus                                                                                                                        |
| 纯化    | Affinity purification                                                                                                                    |

抗原信息


|      |                                                                                                                                                                                                                                                                                                     |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human NAT10.                                                                                                                                                                                                                                            |
| 序列   | KTLTDEDEADQGGWLAAFWKDFRRRFLALLSYQFSTFSPSLALNIIQNRNMGKPAQPALSREELEALFLPYDLKRLEMYSR<br>NMVDYHLIMDMIPAIISRIYFLNQLGDLALSAAQSALLLGIGLQHKSVQDQLEKEIELPSGQLMGLFNRIIRKVVKLFNEVQEKA<br>IEEQMVAAKDVVMEPTMKTLSDDLDEAAKEFQEKHKKEVGKLSMDLSEYIIIRGDDEEWNEVLNKAGPNASIIISLKSDKKR<br>KLEAKQEPKQSKKLKNRETKNKKDKMLKRKK |

靶点信息

|      |                                                                                                                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is an RNA cytidine acetyltransferase involved in histone acetylation, tR<br>NA acetylation, the biosynthesis of 18S rRNA, and the enhancement of nuclear architecture and chromati<br>n organization. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                                  |
|-------|--------------------------------------------------|
| 基因ID  | 55226                                            |
| 基因名   | NAT10                                            |
| Swiss | Q9H0A0                                           |
| 别名    | NAT10; ALP; NET43; Kre33; N-acetyltransferase 10 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)