

# QK1 Rabbit mAb

货号: B30643

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                   |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IF/ICC IP                                                                                  |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50            |
| 理论分子量 | 35kDa/36kDa/37kDa                                                                             |
| 实测分子量 | 38kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HeLa,NIH/3T3,C2C12,C6                                                                         |
| 细胞定位  | Cytoplasm,Nucleus                                                                             |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                                                                                                                                                                                                                                                  |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human QK1.                                                                                                                                                                                                                                                                                                           |
| 序列   | MVGEMETKEKPKPTPDYLMQLMNDKKLMSSLPNFCGIFNHLERLLDEEISVRKDMYNDTLNGSTEKRSAELPDVGPV<br>QLQEKLYVPVKEYPDFNFVGRILGPRGLTAKQLEAETGCKIMVRGKGSMRDKKKEEQNRGKPNWEHLNEDLHVLITVEDA<br>QNRAEIKLKRAVEEVKLLVPAEGEDSLKKMQLMELAILNGTYRDANIKSPALAFSLAATAQAAPRIITGPAPVLPPAALRTP<br>TPAGPTIMPLIRIQITAVMPNGTPHPTAAIVPPGPEAGLIYTPYEYPYTLAPATSILEYPIEPSGVLGAVATKVRRHDMRVHPYQ<br>RIVTADRAATGN |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                              |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is an RNA-binding protein that regulates pre-mRNA splicing, export of mRNAs from the nucleus, protein translation, and mRNA stability. The encoded protein is involved in myel inization and oligodendrocyte differentiation and may play a role in schizophrenia. Multiple transcript vari ants encoding different isoforms have been found for this gene. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                         |
|-------|-------------------------|
| 基因ID  | 9444                    |
| 基因名   | QKI                     |
| Swiss | Q96PU8                  |
| 别名    | QKI;Hqk;QK;QK1;QK3;hqkl |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)