

FOXL2 Rabbit mAb

货号: B30657

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	38kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,Mouse testis,Rat testis,Rat ovary
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human FOXL2.
序列	MMASYPEPEDAAGALLAPETGRTVKEPEGPPSPGKGGGGGGGTAPEKPDPAQKPPYSYVALIAMAIRESAEKRLTSLGIYQYIIAKFPFYEKNKKGWQN

靶点信息

研究背景	This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Expansion of a polyalanine repeat region and o ther mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3. [provided by RefSeq, Jul 2016]
基因ID	668

基因名	FOXL2
Swiss	P58012
别名	BPES; BPES1; PFRK; PINTO; POF3

产品验证



Western blot analysis of FOXL2 expressed in HeLa,K-562,Mouse testis using FOXL2 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn