

CD39 Rabbit mAb

货号: B30664

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa/42kDa/46kDa/57kDa/58kDa/59kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen,Rat lung
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human CD39.
序列	MESEELADRVLDVVERSLSNYPFDQFQGARIITGQEEGAYGWITINYLLGKFSQKTRWFSIVPYETNNQETFGALDLGGAST QVTFVPQNQTIESPDNALQFRLYGKDYNVYTHSFLCYGKDQALWQKLAKDIQVASNEILRDPCHFPGYKKVVNVSDLYKT PCTKRFEMTLPPFQQFEIQG

靶点信息

研究背景	The protein encoded by this gene is a plasma membrane protein that hydrolyzes extracellular ATP and A DP to AMP. Inhibition of this protein's activity may confer anticancer benefits. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	953

基因名	ENTPD1
Swiss	P49961
别名	ENTPD1;ATPDase;CD39;NTPDase-1;SPG64

产品验证



Western blot analysis of CD39 expressed in Mouse spleen,Rat lung using CD39 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn