

# KMT6 / EZH2 Rabbit mAb

货号: B30674

## 产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Monoclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC FC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>FC:</b> 1:20 - 1:50 |
| 理论分子量 | 79kDa/81kDa/84kDa/85kDa/86kDa                                                                                  |
| 实测分子量 | 98kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                  |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | MCF7,293T,HepG2,HeLa,SW620,Raji,Mouse testis                                                                   |
| 细胞定位  | Nucleus                                                                                                        |
| 纯化    | Affinity purification                                                                                          |

## 抗原信息


|      |                                                                                                        |
|------|--------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human KMT6 / EZH2.                                         |
| 序列   | LERTEILNQEWKQRRIQPVHILTSVSSLRGTRECSVTSDLDFPTQVIPLKTLNAVASVPIMYSWSPLQQNFMVEDETVLHNIPYMGDEVLDQDGTFFIEELI |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric p protein complexes, which are involved in maintaining the transcriptional repressive state of genes over successive cell generations. This protein associates with the embryonic ectoderm development protein, the VAV1 oncoprotein, and the X-linked nuclear protein. This protein may play a role in the hematopoietic and central nervous systems. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                                            |
|-------|------------------------------------------------------------|
| 基因ID  | 2146                                                       |
| 基因名   | EZH2                                                       |
| Swiss | Q15910                                                     |
| 别名    | ENX-1;ENX1;EZH1;EZH2b;KMT6;KMT6A;WVS;WVS2;EZH2;KMT6 / EZH2 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)