

KLF4 Rabbit mAb

货号: B30686

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	6kDa/12kDa/45kDa/51kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HCT116,HT-29
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human KLF4.
序列	KRGRRSWPRKRTATHTCDYAGCGKTYTKSSHLKAHLRHTHTGEKPYHCDWDGCGWKFARSDELTRHYRKHTGHRPFQC QKCDRAFSRSDHLALHMKRHF

靶点信息

研究背景	This gene encodes a protein that belongs to the Kruppel family of transcription factors. The encoded zinc finger protein is required for normal development of the barrier function of skin. The encoded protein is thought to control the G1-to-S transition of the cell cycle following DNA damage by mediating the tumor suppressor gene p53. Mice lacking this gene have a normal appearance but lose weight rapidly, and die shortly after birth due to fluid evaporation resulting from compromised epidermal barrier function. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	---

基因ID	9314
基因名	KLF4
Swiss	O43474
别名	EZF;GKLF;KLF4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn