

# CBR1 Rabbit mAb

货号: B30708

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IP                                                                                         |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IP:</b> 1:20 - 1:50                                           |
| 理论分子量 | 18kDa/30kDa                                                                                   |
| 实测分子量 | 30kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | SW480,22Rv1,MCF7,NCI-H460,SKOV3,Mouse kidney,Mouse liver,Mouse lung,Mouse testis              |
| 细胞定位  | Cytoplasm                                                                                     |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                                                                                                                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human CBR1.                                                                                                                                                                                                                                           |
| 序列   | MSSGIHVALVTGGNKGIGLAIVRDLCLRLFSGDVLTARDVTRGQAQVQLQAEGLSPRFHQLDIDDLQSIRALRDFLRKEY<br>GGLDVLVNNAGIAFKVADPTPFHIIQAEVTMKTNFFGTRDVCTELLPLIKPQGRVVNVSSIMSVRALKSCSPELQQKFRSETI<br>TEEELVGLMNKFVEDTKKGVHQKEGWPSSAYGVTKIGVTVLSRIHARKLSEQRKGDKILLNACCPGWVRTDMAGPKATK<br>SPEEGAETPVYLALLPPDAEGPHGQFVSEKRVEQW |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases (SDR) family, w<br>hich function as NADPH-dependent oxidoreductases having wide specificity for carbonyl compounds, suc<br>h as quinones, prostaglandins, and various xenobiotics. Alternatively spliced transcript variants have bee<br>n found for this gene. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                        |
|-------|------------------------|
| 基因ID  | 873                    |
| 基因名   | CBR1                   |
| Swiss | P16152                 |
| 别名    | CBR1;CBR;SDR21C1;hCBR1 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)