

SCF Rabbit mAb

货号: B30712

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	26kDa/27kDa/30kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Rat lung,Rat liver,Mouse testis,Mouse spleen
细胞定位	Cell membrane,Cytoplasm,Secreted,Single-pass type I membrane protein,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SCF.
序列	SKDLKKSFKSPEPRLFTPEEFFRIFNRSIDAFKDFVVASETSDCVSSTLSPEKDSRVSVTKPFMLPPVAA

靶点信息

研究背景	This gene encodes the ligand of the tyrosine-kinase receptor encoded by the KIT locus. This ligand is a pleiotropic factor that acts in utero in germ cell and neural cell development, and hematopoiesis, all believed to reflect a role in cell migration. In adults, it functions pleiotropically, while mostly noted for its continued requirement in hematopoiesis. Two transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4254
基因名	KITLG
Swiss	P21583
别名	KITLG;DCUA;DFNA69;FPH2;FPHH;KL-1;Kitl;MGF;SCF;SF;SHEP7

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn