

Calpain 2 Rabbit mAb

货号: B30718

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	80kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse lung,Mouse brain,Mouse kidney,Rat lung,Rat brain,Rat kidney
细胞定位	Cell membrane,Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Calpain 2.
序列	LWTKIQKYQKIYREIDVDRSGTMNSYEMRKALEEAGFKMPCQLHQVIVARFADDQLIDFDNFVRCLVRLETLFKIFKQLDP ENTGTIELDLISWLCFSVL

靶点信息

研究背景	The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The e mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, r egulatory subunit. This gene encodes the large subunit of the ubiquitous enzyme, calpain 2. Multiple hete rogeneous transcriptional start sites in the 5' UTR have been reported. Two transcript variants encoding d ifferent isoforms have been found for this gene. [provided by RefSeq, Mar 2009]
基因ID	824

基因名	CAPN2
Swiss	P17655
别名	CANP2; CANPL2; CANPml; mCANP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn