

TLE1 Rabbit mAb

货号: B30745

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	83kDa
实测分子量	83kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human TLE1.
序列	DPSSPRASPAHSPRENGIDKNRLKKDASSSPASTASSASSTSLKSKEMSLHEKASTPVLKSSTPTPRSDMPTPGTSATPG LRPGLGKPPAIDPLVNQAA

靶点信息

研究背景	Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. Enhances FOXG1/BF-1- and HES1-mediated transcriptional repression (By similarity. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG.
------	--

基因ID	7088
基因名	TLE1
Swiss	Q04724
别名	TLE1;ESG;ESG1;GRG1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn