

Cyclophilin F Rabbit mAb

货号: B30772

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	16kDa/22kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,B cells,Mouse kidney
细胞定位	Mitochondrion matrix
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Cyclophilin F.
序列	CSKGSGDPSSSSSSGNPLVYLDVDANGKPLGRVVLELKADVVPKTAENFRALCTGEKGFYKYGSTFHRVIPSFMCQAGDFTNHNGTGGKSIYGSRFPDENFTLKHVGPVLSMANAGPNTNGSQFFICTIKTDWLDGKHVVFGHVKEGMDVVKKIESFGSKSGRTSKKIVITDCGQLS

靶点信息

研究背景	The protein encoded by this gene is a member of the peptidyl-prolyl cis-trans isomerase (PPlase) family. P Plases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. This protein is part of the mitochondrial permeability transition pore in the inner m itochondrial membrane. Activation of this pore is thought to be involved in the induction of apoptotic and necrotic cell death.
------	---

基因ID	10105
基因名	PPIF
Swiss	P30405
别名	PPIF;CYP3;CyP-M;Cyp-D;CypD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn