

alpha 1 Catenin Rabbit mAb

货号: B30774

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	59kDa/100kDa/102kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,OVCAR3,293T,MCF7,Mouse lung,Mouse heart,Mouse brain
细胞定位	Cell junction,Cell junction,Cell membrane,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,adhe rens junction,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human alpha 1 Catenin.
序列	ELDDSDFETEDFDVRSRTSVQTEDDQLIAGQSARAIMAQLPQEQKAKIAEQVASFQEEKSKLDAEVSKWDDSGNDIIVLA KQCMIMMEMTDFTRGKGPLKNTSDVISAAKKIAEAGSRMDKLGRTIADHCPDSACKQDLLAYLQRIALYCHQLNICKV KAEVQNLGGELVSGVDSAMSLIQAAKNLMNAVVTQTKASYVASTKYQKSQGMASLNLPVSWKMKAPKPLVKREK QDETQTKIKRASQKKHVNPNVQALSEFKAMDSI

靶点信息

研究背景	This gene encodes a member of the catenin family of proteins that play an important role in cell adhesion process by connecting cadherins located on the plasma membrane to the actin filaments inside the cell. The encoded mechanosensing protein contains three vinculin homology domains and undergoes conformational changes in response to cytoskeletal tension, resulting in the reconfiguration of cadherin-actin filament connections. Certain mutations in this gene cause butterfly-shaped pigment dystrophy.
基因ID	1495
基因名	CTNNA1
Swiss	P35221
别名	CAP102;MDPT2;Catenin alpha 1;CTNNA1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn