

# KLF10 Rabbit mAb

货号: B30777

## 产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                     |
| 理论分子量 | 53kDa                                                                                         |
| 实测分子量 | 53kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HeLa,HepG2,Jurkat,Mouse lung                                                                  |
| 细胞定位  | Nucleus                                                                                       |
| 纯化    | Affinity purification                                                                         |

## 抗原信息


|      |                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human KLF10.                                              |
| 序列   | MLNFGASLQQTAEERMESISRPKESMYSWNKTAEKSDFEAVEALMSMSCSWKSDFKKYVENRPVTPVSDLSEENLLP<br>GTPDFHTIPAFCLTPPYSPD |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of a family of proteins that feature C2H2-type zinc finger domains. The encoded protein is a transcriptional repressor that acts as an effector of transforming growth factor beta signaling. Activity of this protein may inhibit the growth of cancers, particularly pancreatic cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2013] |
| 基因ID | 7071                                                                                                                                                                                                                                                                                                                                                                                                              |

|       |                              |
|-------|------------------------------|
| 基因名   | KLF10                        |
| Swiss | Q13118                       |
| 别名    | EGR-alpha; EGRA; TIEG; TIEG1 |

产品验证



Western blot analysis of KLF10 expressed in Mouse lung,HeLa,HepG2 using KLF10 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)