

# Phospho-IRF3 (S386) Rabbit mAb

货号: B30780

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IF/ICC                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                      |
| 理论分子量 | 12kDa/16kDa/33kDa/47kDa/49kDa                                                                 |
| 实测分子量 | 47kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HeLa,HeLa                                                                                     |
| 细胞定位  | Cytoplasm,Nucleus                                                                             |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                                        |
|------|------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human Phospho-IRF3 (S386). |
| 序列   | ASSLE                                                                  |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the interferon regulatory transcription factor (IRF) family. The encoded protein is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other interferon-induced genes. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. |
| 基因ID | 3661                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

|       |            |
|-------|------------|
| 基因名   | IRF3       |
| Swiss | Q14653     |
| 别名    | IIAE7;IRF3 |

产品验证



Western blot analysis of Phospho-IRF3 (S386) expressed in HeLa,HeLa using Phospho-IRF3 (S386) R abbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)