

Hexokinase 1 Rabbit mAb

货号: B30784

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	101kDa/102kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,MCF7,U-87MG,Mouse testis,Mouse brain,Mouse heart,Rat testis,Rat brain
细胞定位	Mitochondrion outer membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Hexokinase 1.
序列	KIDKYLYAMRLSDETLDIMTRFRKEMKNGLSRDFNPTATVKMLPTFVRSIPDGSEKGDFIALDLGGSSFRILRVQVNHEKN QNVHMESEVYDTPENIVHGSQSLFDHVAECLGDFMEKRKIKDKKLPVGFTFSFPCQSQSKIDEAILITWTKRFKASGVEGA DVVKLLNKAIKRGDYDANIVAVVNDTVGTMTCGY

靶点信息

研究背景	Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in this gene have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of this gene results in several transcript variants which encode different isoforms, some of which are tissue-specific.
------	---

基因ID	3098
基因名	HK1
Swiss	P19367
别名	HK1;HK;HK1-ta;HK1-tb;HK1-tc;HKD;HKI;HMSNR;HXX1;hexokinase;RP79

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn