

Prohibitin Rabbit mAb

货号: B30790

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	17kDa/29kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,293T,Neuro-2a,Mouse liver,PC-12,Rat liver
细胞定位	Mitochondrion inner membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Prohibitin.
序列	MAAKVFESIGKFGGLALAVAGGVVNSALYNVDAGHRAVIFDRFRGVQDIVVGEGTHFLIPWVQKPIIFDCRSRPRNVPVITGS KDLQNVNITLRLFRPVASQLPRIFTSIGEDYDERVLPSITTEILKSVVARFDAGELITQRELVSRQVSDDLTERAATFGLILDD VSLTHLTFGKEFTEAVEAKQVAQQAERARFVVEKAEQKKAIIISAEGDSKAAELIANSLATAGDGLIELRKLEAAEDIAYQ LSRSRNITYLPAGQSVLLQLPQ

靶点信息

研究背景	This gene is evolutionarily conserved, and its product is proposed to play a role in human cellular senescence and tumor suppression. Antiproliferative activity is reported to be localized to the 3' UTR, which is proposed to function as a trans-acting regulatory RNA. Several pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants.
基因ID	5245
基因名	PHB
Swiss	P35232
别名	HEL-215;HEL-S-54e;PHB1;PHB;Prohibitin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn