

ADRA1A Rabbit mAb

货号: B30798

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	32-52kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,LO2,HL-60,Mouse heart,Mouse kidney,Mouse spinal cord,Rat heart
细胞定位	Cell membrane,Multi-pass membrane protein,Nucleus membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ADRA1A.
序列	LCRKQSSKHALGYTLHPPSQAVEGQHKDMVRIPVGSRETFYRISKTDGVCEWKFFSSMPRGSARITVSKDQSSCTTARTK SRSVTRLECSGMILAHCNLRLPGSRDSPASASQAAGTTGMCHQADATRPS

靶点信息

研究背景	Alpha-1-adrenergic receptors (alpha-1-ARs) are members of the G protein-coupled receptor superfamily. They activate mitogenic responses and regulate growth and proliferation of many cells. There are 3 alpha-1-AR subtypes: alpha-1A, -1B and -1D, all of which signal through the Gq/11 family of G-proteins and different subtypes show different patterns of activation. This gene encodes alpha-1A-adrenergic receptor. Alternative splicing of this gene generates four transcript variants, which encode four different isoforms with distinct C-termini but having similar ligand binding properties.
------	--

基因ID	148
基因名	ADRA1A
Swiss	P35348
别名	ADRA1A;ADRA1C;ADRA1L1;ALPHA1AAR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn