

Perilipin A Rabbit mAb

货号: B30804

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	55kDa
实测分子量	62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	cytosol,endoplasmic reticulum,lipid droplet
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Perilipin A.
序列	HYVPLPRLSLMEPESEFRDIDNPPEAVERREAERRASGAPSAGPEPAPRLAQPRRSLRSAQSPGAPPGGLEDEVATPAAPRPGFPAVPREKPKRRVSDSF

靶点信息

研究背景	The protein encoded by this gene coats lipid storage droplets in adipocytes, thereby protecting them until they can be broken down by hormone-sensitive lipase. The encoded protein is the major cAMP-dependent protein kinase substrate in adipocytes and, when unphosphorylated, may play a role in the inhibition of lipolysis. Alternatively spliced transcript variants varying in the 5' UTR, but encoding the same protein, have been found for this gene.
基因ID	5346

基因名	PLIN1
Swiss	O60240
别名	PLIN1;FPLD4;PER1;PLIN

产品验证



Western blot analysis of Perilipin A expressed in HepG2 using Perilipin A Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn